A bed grab rail is designed to give assistance and support:

- When raising yourself from a lying to sitting position and as you stand from the bed.
- Whilst you are stood beside the bed, for example whilst getting dressed.
- As you change your position in bed.
- As you sit and lift your legs up onto the bed.



For qualified, free*
professional advice
and support,
Sue and her team
are ready to help.

*includes a big smile





Orderline

0845 260 7061

Buy online

www.otstores.co.uk

Email info@otstores.co.uk


OTS Ltd The Haysfield, Spring Lane North, Malvern WR14 1GF

bed grab rails



Let a bed grab rail take the strain when you get up

Choose the best one for you



This increasingly popular all-metal bed grab rail is simple to assemble and dismantle without tools, making it an ideal choice if you want to take it with you when going on holiday, or visiting relatives. Flat when dismantled for easy storage. Easy to clean. Designed for use on divan style beds; suitable for use on single, double or kingsize beds.

An additional advantage of the Easyleaver® range is the choice of sizes depending on the user's weight. The standard Easyleaver® bed grab rail has a maximum user weight of 17 stone (108kg), with a Heavy Duty version accommodating a maximum user weight of 25 stone (160kg).

Choice of two colours – Classic White or Silver-grey.

For those over 25 stone ask about our new Bariatric version.

The Liftwell Bed Leaver® bed grab rail is probably the most well recognised bed grab rail on the market, used widely by health professionals across the NHS and Social Services. The baseboard is made from high grade ply, which will not crack or split, and the beech hardwood end block, which is visible on the outside edge of the bed when fitted, provides an attractive finish.

Easyrail® Bed Grab Rail - Single Handle

Easyrail® single-handle can be used on single, double or kingsize beds, and has one handle.

Choice of sizes determined by user's weight – standard Easyrail® has a maximum user weight of 20 stone (127kg), whilst the Heavy Duty version has a maximum user weight of 30 stone (190kg).

For those over 30 stone ask about our new Bariatric version.

Specially designed to fit those beds where the mattress sits within a side frame, such as pine slatted beds.

Easy to assemble and dismantle, with no need for tools.

Can easily be transported or stored when not in use.

Choice of 2 colours for both Standard and Heavy Duty sizes – classic White, or new Silver-Grey.

Easyrail® Bed Grab Rail, Twin-Handled

If two occupants need assistance, or a single occupant would benefit from the use of a bed grab rail handle on both sides of the bed, the Easyrail® twin-handled is also available.



Choosing an OTS bed grab rail

Getting in and out of bed can be quite a struggle for many of us. Sometimes a little support is all that we need to get up and going.

The chart below helps you to determine the most appropriate bed rail in relation to the type of bed you have, and the user's weight. Please contact advice@otstores.co.uk or 0845 260 7061 if you need any further advice.

Type of bed ▶	Divan style bed (Single, double or kingsize)	Side-framed / slatted bed (Single, double or kingsize)
Weight of user ▼ Up to 17 Stone	Easyleaver® Standard (see EASY217 or EASY218) or Bed Leaver® (see BEDL200)	
Up to 20 Stone		Easyrail® Standard (see RAIL250 or RAIL251) Available with single or twin handles
Up to 25 Stone	Easyleaver® Heavy Duty (see EASY225 or EASY226)	
Up to 30 Stone		Easyrail® Heavy Duty (see RAIL350 or RAIL351)
Up to 50 Stone	Easyleaver® Bariatric (see EASY290)	Easyrail® Bariatric (see RAIL390)



Orderline

0845 260 7061

Buy online

www.otstores.co.uk

Email info@otstores.co.uk